IVBPA

TENSIÓMETRO DE BRAZO INALÁMBRICO
WIRELESS ARM BLOOD PRESSURE MONITOR
TENSIÔMETRO DE BRAÇO SEM FIO
TENSiomètre pour bras sans fils

The smartphone screen shows the app's main interface with a user profile picture and name. Below it, a list of recent measurements is displayed:

Medida	Data	Sistólico	<th freq.<="" th=""></th>	
BÁSCULA	22 Jul 2015	64.3	22.1	PPM
BRAZO	17 Jul 2015	63/113	MMHG	62 PPM
MUÑECA	19 Jul 2015	60/111	MMHG	59 PPM

On the right, a table of contents lists the manual versions:

ES	GUÍA DE USUARIO	2
EN	USER GUIDE	10
PT	MANUAL DE INSTRUÇÕES	17
FR	GUIDE DE L'UTILISATEUR.....	24

The blood pressure monitor device is shown at the bottom, with its digital display showing systolic, diastolic, and pulse rates. The device has a white strap and a control button labeled "START/STOP".

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Precision: +-3 mmHg/g

Sensitivity: 1 mmHg

Blood pressure measurement range: 0 - 280 mmHg

Pulse rate measurement range: 40 - 199 beats/min

Temperature/ humidity for use:

+5° C to +40° C / max: 15 to 80 % HR

Temperature/ humidity for storage:

-20° C to +50° C / max: 15 to 90 % HR

On-line and off-line operation

INTENDED USE

The IVY HEALTH IVBPA blood pressure device allows you or your doctor to take your blood pressure using non-invasive means. All you have to do to charge the device's internal battery is connect the cable included to a charging socket or PC USB port. It's that simple: the blood pressure device is then ready to take your blood pressure and pulse rate.

SAFETY WARNINGS

Do not use this device for self-medication on the basis of its measurements or for any therapeutic measure without consulting your doctor.

Contraindications:

1. This blood pressure monitor should not be used if you:
 - a. Suffer from severe arrhythmia, pre-eclampsia, ventricular and/or arterial occlusive disease
 - b. Are undergoing intravascular therapy
 - c. Have venous access
 - d. Have had a mastectomy
 - e. Are pregnant
2. Do not allow its use by children.
3. If you suffering from any illness, consult your doctor before using this blood pressure monitor.

Safety warnings on the battery and other parts of the device:


1. Never open the blood pressure device or dismantle the batteries.
2. Keep the batteries out of reach of children. Do not throw them in fire or short circuit them as they could explode.

3. Do not throw the batteries or any of the device components in ordinary rubbish. Take them to specialised collection centres or put them in a proper container if there is one in your area.
4. Do not use the blood pressure device if you are allergic to any of its components.
5. Do not use the blood pressure device for purposes other than those indicated in this manual

Cleaning the device:

1. Use a damp cloth with a soap solution, but make sure no liquid gets into the device. The blood pressure device should not be used if it is wet.
2. Do not use abrasive liquids or cloths or hard utensils when cleaning.
3. It's recommended to disinfect the inside of the sleeve in contact with the skin with a soft cloth slightly moistened with ethyl alcohol from time to time.
4. Do not inflate it until it has been correctly positioned on your arm.
5. Do not expose it directly to the sun, extreme humidity and/or dust.

LED INFORMATION ON SCREEN

- ① Diastolic blood pressure
- ② Systolic blood pressure
- ③ Pulse rate per minute
- ④ Heartbeat detection

- ⑤ "START/STOP" button
- ⑥ Battery indicator
- ⑦ Bluetooth signal indicator

Some of the possible error messages:

Err 1 Two measurements were performed in a very short space of time and the sleeve has not had enough time to deflate properly. Remove the arm strap, wait a few seconds and try to repeat the measurement.

Err 2 The measurement was not successful. Please try again.

Err 4 The sleeve is not adjusted properly and the electronic inflation cannot be executed. Please place again the armband correctly as explained below and get the measurement repeated.

Err 6/7 Error detecting diastolic pressure. The cuff was moved during the measurement process. Please try again.

Err0r Bluetooth communication generic error.

Err0r LO Low battery. Please charge the device.


WHAT CAN I DO WITH THIS BLOOD PRESSURE DEVICE?

Can it work without the application and my mobile phone?

Yes, you can use the blood pressure device without the free IvyHealth app and without connection to a mobile phone or tablet. You can see the results on the device's screen. However, if you are connected to a mobile phone with the IVY Cloud app open, you will have a more complete view of your blood pressure and other values.

How do I get my blood pressure measured?

The best way is to sit comfortably. Put your arm at a right angle, place the cuff on your arm with the red "UP" arrow facing up. The bottom edge of the cuff should be 1 to 2 cm above the elbow. You are now ready to go ! Just press the "START / STOP" button and wait for a few seconds to get your first measurement.

You are recommended to take the readings at the same time and under the same conditions in order to be able to compare results. And it is always better to take your blood pressure before meals and after resting for a few minutes.

What if I download the application?

Download the free IvyHealth application from the Apple Store or Google Play store. When you first open the application, you must register and set up your user name. It is very easy if you follow the instructions on your mobile phone's screen.

The IvyHealth application can be automatically paired with your recently acquired blood pressure device.

1. When you first open the application you'll find the registration screen. Once you are registered you will enter the main menu, where you'll be able to choose the device you'll be working with.
2. To make a measurement, press the icon on the main screen. You will access the measurement screen and can now position the blood pressure monitor on your arm and follow the instructions on the screen.
3. Press the START/STOP key on your blood pressure device.
4. The screen of the device will show the data corresponding to the provisional reading. Your phone will automatically detect the measurement and will show the collected data.
5. Wait until measurement has been completed.
6. Once the transmission has been completed, you can see the following data (*) on your mobile phone screen:
 - a. Your systolic pressure
 - b. Your diastolic pressure
 - c. Your pulse rate
7. Once the measurement is finished, click on the "report" button and you will see a highly intuitive report that tells you if the readings made are at the required level or if your blood pressure is too high or too low, in accordance with the readings laid down by the WHO.
7. Press to save and your first reading is recorded for the future. You can consult your phone or tablet whenever you wish. You will be protected by your own password.

USER GUIDE

8. Subsequent measurements will be kept in your mobile phone whenever you measure your blood pressure with the app open and you follow the previous steps.

What else can I do with the application?

If you use the blood pressure monitor together with the free IVY Cloud application, you can obtain the readings for all the body values listed in the previous paragraph.

However, the application also allows you to:

1. Check your tendencies. Just click on the  icon, and then access to "history".
2. Save and display the results of as many users as you wish. Every time somebody makes a measurement from their mobile phone, and if they are registered as a user, they can record their results and protect them with a password

In addition to blood pressure

The IVY HEALTH package of health complements is made up of different devices affording you greater control over your health.

Check out our latest healthcare products following this link:
<http://www.ivyhealth.es>

LEGAL NOTICE

Declaration of Conformity

FONEXION SPAIN S.A. declares that the arm blood pressure monitor IVBPA is on compliance with the european regulations EN 1060-1 and EN 1060-4, in addition with the essential requirements and other relevant provisions of Directives MDD93/42/CE and R&TTE 1999/5/CE.

You can get the full declaration of conformity here:

<http://www.fonexion.com/doc>.

Warranty

This product has the warranty period which is legally applicable in each country, starting from the date of purchase. Save your sales receipt, is your warranty certificate. Please keep the original packaging and components to execute the warranty, handling it on the place of purchase. The warranty shall not be valid on the following cases: documents which are illegible or voluntarily modified, where the product has been misused or modified by unauthorised personal; or damaged owing to misuse or negligence, in addition to any damage deriving from storage or shipment thereof.

Bluetooth

The Bluetooth word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by IVY HEALTH is under license. Any other trademarks and/or trade names are those of their respective owners.

 **STANDARDS AND DOCUMENTS**

Title	Code
Medical Device Quality System for the regulation requirement	EN ISO 13485:2012
Medical device risk management to medical devices application	EN ISO14971:2012
Symbol for the label of medical devices	EN 980:2008
Term, symbol and information of medical device— information of medical device manufacturer offering	EN 1041:2008
Medical electrical equipment – Part 1: General requirements for safety	EN 60601-1: 2005/AC:2012
Medical electrical equipment – Part 1-2: General requirements for safety – Collateral standard: Electromagnetic compatibility – Requirements and tests	EN 60601-1-2:2007
Clinical investigation of medical devices for human subjects. General requirements	EN ISO 14155-1:2009
Clinical investigation of medical devices for human subjects. Clinical investigation plans	EN ISO 14155-2:2009
Medical device software – Software life-cycle processes	EN 62304:2006/AC: 2008
Medical electrical equipment – Part 1-6: General requirements for basic safety and essential performance – Collateral standard: Usability	EN 60601-1-6:2010
Medical electrical equipment. Part 1-11:General requirements for basic safety and essential performance. Collateral standard:Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment.	EN 60601-1-11:2010

USER GUIDE

Medical devices – Application of usability engineering to medical devices	EN 62366: 2008
Council Directive 93/42/EEC of 14 June 1993 concerning medical devices OJ L 169 of 12 July 1993 (including Directive 2007/47/EC)	MDD 93/42/EEC
Non-invasive sphygmomanometers —Part 1: General requirements	EN 1060-1
Non-invasive sphygmomanometers. Part 3:Supplementary requirements for electro-mechanical blood pressure measuring systems	EN 1060-3
Secondary cells and batteries containing alkaline or other non-acid electrolytes – Safety requirements for portable sealed secondary cells, and for batteries made from them, for use in portable applications	IEC 62133:2012
User Guide	Versión 1.0

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FONEXION SPAIN S.A. GRAN VÍA 45, 1º PLANTA DPT.114. 48011 BILBAO
Manufactured for **FONEXION SPAIN S.A.**
EC REP Wellkang Ltd · Suite B, 29 Harley Street · London, W1G 9QR, U.K.
Shenzhen Belter Health Measurement and Analysis Technology CO., Ltd
Fabricado para FONEXION SPAIN S.A.